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## LIST OF CONTENTS

Items Supplied	9
Controller – At A Glance	9
Symbols	10
Signal Words	10
Intended Use	10
When May the Device Not Be Used?	11
Safety Notices	11
Inserting / Replacing Batteries	12
Safety Devices	12
Overview of Functions and Programmes	12
Use	13
Cleaning and Storage	14
Troubleshooting	15
Disposal	15
Technical Data	16

## ITEMS SUPPLIED

(Picture A)

- 1 Thigh / buttock pad 1 x
  - 2 Controller 4 x
  - 3 Arm / leg pad 2 x
  - 4 Abdomen pad 1 x
- Operating instructions 1 x

Check the items supplied for completeness and the components for transport damage. If you find any damage, do not use the device but contact our customer service department.

Remove any possible films, stickers or transport protection from the device.

**Never remove the rating plate and any possible warnings!**

Two batteries (1.5 V; AAA) per controller (2) are also required.

## CONTROLLER – AT A GLANCE

(Picture B)

- 5 Programme lamps
- 6 OFF / – button
- 7 MODE button
- 8 ON / + button

**Dear Customer,**

We are delighted that you have chosen the MAXXMEE EMS set.

The comfortable and yet effective method of EMS (electromagnetic muscle stimulation) produces targeted muscle contractions. Thanks to the self-adhesive gel film, the pads are simply stuck onto the particular region of the body and enable targeted stimulation of the muscles.

If you have any questions about the device and about spare parts/accessories, contact the customer service department via our website: [www.dspro.de/kundenservice](http://www.dspro.de/kundenservice)

We hope you have a lot of fun with your MAXXMEE EMS set.

**Information About the Operating Instructions**

Before using the device for the first time, please read through these operating instructions carefully and keep them for future reference and other users. They form an integral part of the device. The manufacturer and importer do not accept any liability if the information in these operating instructions is not complied with.

**SYMBOLS**

Danger symbol: This symbol indicates possible dangers. Read the associated safety notices carefully and follow them.

Supplementary information

Read operating instructions before use!

Symbol for direct current

**SIGNAL WORDS**

- WARNING** warns of **possible** serious injuries and danger to life
- NOTICE** warns of material damage

**INTENDED USE**

- The device is intended to be used to stimulate the muscles of the abdomen, legs, arms, back and buttocks.
- Do **not** use the device while driving a vehicle or while sleeping.
- The device should **not** be used before going to sleep. The stimulation may cause you problems getting to sleep.
- The device is not designed for medical or therapeutic purposes.
- The device is for personal use only and is not intended for commercial applications.
- The device should only be used for the specified purpose and only as described in the operating instructions. Any other use is deemed to be improper.
- The warranty does not cover faults caused by incorrect handling, damage or attempted repairs. The same applies to normal wear and tear.

## WHEN MAY THE DEVICE NOT BE USED?

- This device is not intended to be used by people (including children) with reduced sensory or mental abilities or a lack of experience and/or knowledge, unless they are supervised by a person who is responsible for their safety or have received instruction from this person on how to use the device safely. There is a danger of injury from incorrect use!
- Infants and people who are sleeping or unconscious must not be treated with the device.
- The following people may use the device only following consultation with a doctor:
  - people who are fitted with a pacemaker or similar medical implants;
  - people with heart disease or who suffer from epilepsy;
  - women who are pregnant or have just given birth;
  - people with varicose veins;
  - people suffering from diabetes;
  - people who have a malignant tumour;
  - people who suffer from circulation problems;
  - people with nerve damage or paralysis;
  - people who belong to a risk group or who have concerns about their health.
- If you are unsure whether the device may be used, consult a responsible doctor.
- Do not use the device on:
  - body parts with no feeling;
  - regions which have recently been operated on or are scarred;
  - the immediate vicinity of injured, infected, inflamed, (sun-)burnt or very sensitive skin;

- swelling;
- skin irritations, acne, birthmarks, warts, rosacea;
- a skin rash or fungal infections.
- Do not use the device on your head, neck, face, genital area or on the rib cage close to your heart.
- Do not use the device if you have calf pain.
- Do not use the device if you have a fever.

## SAFETY NOTICES

### WARNING – Danger of Injury

- Keep children and animals away from plastic bags and films. There is a danger of suffocation.
- Children must be supervised to ensure that they do not play with the device. There is a danger of injury from incorrect use!
- Do not use the device in the immediate vicinity of devices which operate using short waves or microwaves. Maintain a minimum gap of 1 metre.
- The device must not be used near a mobile phone. Do not use a mobile phone while using the device.
- Do not use the device if you are attached to high-frequency surgical devices or are in the vicinity of such devices or high-frequency commercial equipment.
- Batteries can be deadly if swallowed. Batteries and controllers should be kept out of the reach of children and animals. If a battery is swallowed, seek medical assistance immediately. There is a danger of burns from battery acid!
- If a battery leaks, do not allow your skin, eyes and mucous membranes to come into contact with the battery acid. If you touch battery acid, rinse the affected areas with abundant clear

water straight away and seek medical assistance immediately. There is a danger of burns from battery acid!

- Do not use the device if it is damaged.

### NOTICE – Risk of Damage to Material and Property

- Use only the battery type specified in the technical data.
- Always replace all batteries together. Never use batteries of different types or brands, or batteries with different capacities.
- Take the batteries out of the controllers when they are dead or if you will not be using the device for a prolonged period. This will prevent any damage that might occur if a battery leaks.
- Remove a dead battery from the controller immediately.
- Do not use the device in rooms with high humidity.
- Do not bend the pads sharply.
- Switch off the controllers before changing the position of the pads.
- The pads must not come into contact during use with metal parts, e.g. on belts, wrist watches, jewellery or similar items.
- Protect the device from: impacts, direct sunlight, sub-zero temperatures, persistent moisture and wet conditions.
- Use only original accessories from the manufacturer in order to guarantee that there is no interference that may prevent the device from working and to avoid possible damage.

## INSERTING / REPLACING BATTERIES

Each controller (2) requires 2 batteries (1.5 V; AAA).

Before inserting the batteries, check whether the contacts in the controller and on the batteries are clean, and clean them if necessary.

1. Slide open the battery compartment cover on the bottom of a controller in the direction of the arrow.
2. Remove any used batteries and dispose of them in an environmentally friendly way.
3. Insert two new batteries. Pay attention to the correct polarity (+ / -).
4. Place on the battery compartment cover and slide it counter to the direction of the arrow. The locking lug must engage audibly.

## SAFETY DEVICES

The device has two built-in safety functions:

- The device automatically switches off after 10 seconds if no contact is established between the skin and silicon electrode.
- The device automatically switches off after 15 minutes of use.

## OVERVIEW OF FUNCTIONS AND PROGRAMMES

### Functions

- The ON / + button (8) is used to switch on the device, set the intensity of the muscle stimulation and confirm the selected programme. The controllers have 9 intensity levels. The controllers start on level 1. If you keep pressing the button, the next-higher

intensity level is set. The button needs to be pressed 9 times to reach the maximum intensity level.

- The OFF / – button (6) is used to switch off the device or decrease the intensity level. Press the button once to set the lower intensity level in each case.
- The MODE button (7) is used to select a stimulation programme. The programme lamps (5) indicate which programme has been selected (see "Programmes").

### Intensity Level

Level	Intensity	Function
1 – 3	Weak	Warming up
4 – 7	Moderate	Training
8 – 9	Strong	Fat burning

### Programmes

You can press the MODE button (7) to switch between the following stimulation programmes:

Programme/ Lamps	Function	Description
1 ●○○○	Gentle warming up	Alternate between 3 seconds of muscle contractions and 2 seconds of a break
2 ●●○○	Warming up / warming down	Rapid muscle contractions, tapping massage
3 ○○●○	"Running" training mode	Stimulation with fast and slow elements

Programme/ Lamps	Function	Description
4 ○○○●	"Swimming" training mode	Alternating routine with lots of fast pulses
5 ●●○○	"Increasing consumption" training mode	Very many, rapid muscle contractions
6 ○○●●	"Muscle building" training mode	Constant stimulation with long muscle contractions

### USE

#### Please Note!

- For reasons of hygiene, a pad should only be used by one person in order to prevent the transfer of skin infections or similar complaints.
- If you experience an unpleasant tingling sensation while using the device, stop using it immediately!
- If you experience increased pain or increased tension, immediately switch off the device and stop using it and consult a doctor.
- Do not use straight after eating.
- We explicitly point out that with some people skin irritations or hypersensitivity may occur as a result of the electrical muscle stimulation. In this case, stop using the device immediately and consult a doctor if necessary.

Beginners should stimulate a muscle group for no more than 15 minutes per day! We recommend beginners do no more than 1 – 2 units per day with a rest phase of at least 4 hours between the units!

Advanced users should stimulate each muscle group for no more than 30 minutes a day. We recommend no more than 3 units per day, with a rest phase of at least 4 hours between the units. We explicitly point out that improper use of the device and excessive muscle stimulation can be detrimental to your health.

1. Attach a controller (2) to a pad using the press studs.
2. Clean the areas of skin that are to be stimulated beforehand.
3. Peel the protective film off the pad and place the silicon electrodes on the body so that the appropriate muscle groups can be stimulated (e.g. shoulders, upper arm, lower arm, back, abdomen, hips, thighs or calves). The optimum position is where the strongest contractions are felt.

 Tip: You can stimulate several muscle groups simultaneously by positioning all pads at different places.

4. Once a pad has been placed on the muscle group that is to be stimulated, switch on the controller using the ON / + button (8).
5. Use the **Mode** button (7) to select a programme (see the "Overview of Functions and Programmes" chapter). The corresponding programme lamp(s) (5) light up.
6. Confirm the chosen programme with the ON / + button. The corresponding programme lamp(s) flash.

 If the programme lamp(s) only flash briefly and then light up constantly, the pads are not properly in contact with the skin and the device switches off after approx. 10 seconds.

7. Use the ON / + and OFF / – buttons (6) to set the intensity level you want (1 – 9) (see the "Overview of Functions and Programmes" chapter).

 The device automatically switches off after 15 minutes. To switch it off before this, press the OFF / – button several times if necessary.

8. Take off the pad and clean it after each use in the manner described in the "Cleaning and Storage" chapter.

## CLEANING AND STORAGE

### Please Note!

- No liquid must get into the controllers. Never immerse the controllers in water!
- Do not use any sharp or abrasive cleaning agents or cleaning pads to clean the device. These may damage the surfaces.

1. Remove the controllers (2) from the pads after use.
2. Clean the used pads under running water and allow them to dry in the air.  
**Do not use soap!**
3. Clean the controllers with a slightly damp cloth if necessary.
4. Cover the silicon electrodes with the protective films when the pads are completely dry.
5. Store all parts in a dry place which is protected from moisture, dirt and direct sunlight and is out of the reach of children and animals.

## TROUBLESHOOTING

If the device does not work properly, check whether you are able to rectify the problem yourself. If the problem cannot be solved with the steps below, contact the customer service department.

**Do not attempt to repair an electrical device yourself!**

**A controller (2) does not work. / The programme lamps (5) do not light up.**

No batteries or empty batteries have been inserted in the controller.

- ▶ Insert / replace batteries.

The batteries are inserted incorrectly.

- ▶ Pay attention to the correct polarity when inserting the batteries.

**The programme lamps (5) light up, but do not flash.**

There is no contact between the skin and silicon electrodes.

- ▶ Make sure that the pads fit correctly.

The controller (2) is not correctly attached to the pad.

- ▶ Check whether the controller is correctly attached to the press studs.

**The programme lights (5) are flashing, but no pulses are being sent to stimulate the muscles.**

The pad is not placed on correctly.

- ▶ Keep changing the position of the pad until the spot with the strongest muscle contraction has been found.

The intensity level is set too low.

- ▶ Press the **ON / +** button (8) to increase the intensity level.

**The pads no longer stick.**

The pads are not clean.

- ▶ Clean the used pads under running water.

The silicon electrodes are too old.

- ▶ A replacement set can be ordered via our customer service department.

## DISPOSAL

Dispose of the packaging material in an environmentally friendly manner so that it can be recycled.

Batteries should be removed from the controllers before disposing of the device and disposed of separately from the device at suitable collection points.

This device is governed by the European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). Do not dispose of this device as normal domestic waste, but rather in an environmentally friendly manner via an officially approved waste disposal company.

## TECHNICAL DATA

Article number:	08217
Model number:	TV-W191128
Power supply:	each controller 3 V (2 x 1.5 V battery; AAA)
Frequency range:	1 – 100 Hz
Operating time:	15 minutes
ID of operating instructions:	Z 08217 M DS V1 0320 uh

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